

## Assay of citrate synthase

### Principle:

- Citrate synthase is an enzyme of the tricarboxylic acid cycle (Krebs cycle). It catalyses the formation of citrate from oxaloacetate and acetylCoA. Citrate synthase activity is often used to evaluate the mitochondrial mass in a given tissue. The reduced CoA (CoA-SH) formed during the reaction transforms the 5,5' dithiobis 2 nitrobenzoic acid (DTNB) into 2-nitro-5-benzoic acid (TNB), which absorbs specifically at 412 nm. Citrate synthase activity will therefore be assessed by following the increase of TNB absorbance at 412 nm.
- A baseline activity exists and must be subtracted from the specific activity.

### Practical set up:

- 6 samples may be analyzed on the same day (5 patients and 1 control).
- Reaction medium composition:
  - 100  $\mu$ M DTNB
  - 100 mM Tris HCl pH 8.1
  - 300  $\mu$ M acetylCoA
  - 500  $\mu$ M oxaloacetate
  - 0.1 % Triton X100
  - Tissue: 40  $\mu$ g of proteins (post-nuclear supernatant or cell suspension) or 4  $\mu$ g of mitochondrial proteins
- Preparation of reaction medium:
  - 1) In a 15 mL tube, prepare enough reaction medium for 7 samples i.e. 7 assays:

Reagents	Global quantity	Quantity/sample
5 mM DTNB	140 $\mu$ L	20 $\mu$ L
10 mM acetylCoA	210 $\mu$ L	30 $\mu$ L
10% Triton X-100	70 $\mu$ L	10 $\mu$ L
1 M Tris HCl pH 8.1	665 $\mu$ L	95 $\mu$ L
H <sub>2</sub> O	5425 $\mu$ L	775 $\mu$ L

- 2) ) In a 1 mL cuvette, add **930  $\mu$ L of the reaction medium and 20  $\mu$ L of post nuclear supernatant or skin fibroblasts (diluted at a final protein concentration of 2 mg/mL) or of isolated beef heart mitochondria diluted at 1/200.** Mix.
- Assay:
    - 1) Reading in the spectrophotometer, at 37°C, at wavelength 412 nm  
Initial calibration is performed on air.
    - 2) Incubate the cuvettes **at 37°C, during 5 min**, in the spectrophotometer
    - 3) Read baseline every 20 seconds during 4 minutes  
Six cuvettes may be read at the same time.
    - 4) Start the reaction with **50  $\mu$ L of 10 mM oxaloacetate in 100 mM Tris HCl pH 8.1**
    - 5) Read the reaction every 20 seconds during 4 minutes  
Six cuvettes may be read at the same time.  
If the increase of absorbance appears too rapid giving a non-linear curve, re-do the assay with less tissue.
  - Calculation:
    - 1) Citrate synthase activity is calculated as nanomoles/min/mg proteins.
    - 2) The extinction coefficient for DTNB is  $\epsilon = 13.6$
    - 3) The correction factor is therefore 1838.2 for 40  $\mu$ g of proteins and 18382 for 4  $\mu$ g of proteins in the assay.