

## Assay of Complex II (succinate ubiquinone oxido-reductase)

### Principle:

- Complex II in the respiratory chain transfers electrons from succinate to ubiquinone. Its activity will be assessed by following the reduction of 2,6-dichlorophenolindophenol (DCPIP) by the decrease of the absorbance at 600 nm of the oxidized DCPIP. Ubiquinone (or coenzyme Q<sub>10</sub>), the very hydrophobic natural acceptor, is replaced by decylubiquinone, a more hydrophilic component.
- Baseline is subtracted from the activity started by the addition of decylubiquinone to take into account any non-specific electron transfer to DCPIP.
- The specific complex II activity is sensitive to malonate. Sensitivity to malonate is usually above 90% of the activity observed without malonate. Activity in the presence of malonate is therefore not useful in routine diagnostic procedure.

### Practical set up:

- 6 samples may be analyzed on the same day (5 patients and 1 control). This means 6 assays.
- Reaction medium composition:
  - 25 mM K Phosphate pH 7.5
  - 20 mM succinate
  - 100 µM decylubiquinone
  - 50 µM DCPIP
  - 1 mM KCN
  - 100 µM ATP
  - 2 mg/mL BSA
  - Tissue: 40 µg of proteins (post-nuclear supernatant (liver or muscle) or cell suspension) or 4 µg of mitochondrial proteins
  - Inhibition may be performed with malonate but is not mandatory for diagnosis.

- Preparation of reaction medium:

1) In a 15 mL tube, prepare enough reaction medium for 7 samples i.e. 7 assays:

Reagents	Global quantity	Quantity/sample
500 mM K Phosphate pH 7,5	350 µL	50 µL
50 mg/mL BSA	280 µL	40 µL
200 mM succinate	700 µL	100 µL
5 mM DCPIP	70 µL	10 µL
10mM KCN	700 µL	100 µL
10mM ATP	70 µL	10 µL
H <sub>2</sub> O	4662 µL	666 µL

2) In a 1 mL cuvette, add **976 µL of the reaction medium and 20 µL of post nuclear supernatant or cell suspension** (diluted at a final protein concentration of 2 mg/mL) **or of isolated beef heart mitochondria** (diluted at 0.2 mg/mL i.e. 1/200). Mix.

- Assay:

- 1) Reading in the spectrophotometer, at 37°C, at wavelength 600 nm  
Initial calibration is performed on air.
- 2) Incubate the cuvettes **at 37°C, during 5 min**, in the spectrophotometer
- 3) Read baseline every 15 seconds during 3 minutes
- 3) Start the reaction by adding **4 µL of 25 mM decylubiquinone kept at room temperature**.
- 4) Reading every 15 seconds during 3 minutes,  
If the decrease appears too rapid giving a non-linear curve, re-do the assay with less tissue.

- Calculation:

- 1) Complex II specific activity is calculated as nanomoles/min/mg proteins.
- 2) The extinction coefficient for DCPIP is  $\epsilon = 19.2$
- 3) The correction factor is therefore 1302.1 for 40 µg of proteins and 13021 for 4 µg of proteins in the assay.