

Tissue preparation for the spectrophotometric assays of the respiratory chain activities

Principle:

- All samples (muscle, liver, peripheral blood mononuclear cells, cultured skin fibroblasts) have to have been quickly frozen and continuously kept at -80°C afterwards. Shipment has to be done in dry ice.
- Tissue preparation is performed at $0-6^{\circ}\text{C}$.

Protocol:

- 5 patients samples may be analyzed in one day. An aliquot of the beef heart mitochondria is always analyzed in parallel to serve as an internal control of the assay quality.
- **2 series** of 5 1.5 mL Eppendorf tubes are numbered and kept in ice
- **Preparation of 5% homogenates in "Mannitol buffer" (muscle, liver fragments):**
 - 1) 5 glass Potter homogenizers (Kontes, size 20) are numbered and kept in ice.
 - 2) 20 to 50 mg tissue fragment are used
If needed, cut the fragment without unfreezing it:
 - either with a scalpel, the fragment being held with tweezers in a plastic box, both tweezers and plastic box cooled to -80°C by their contact with dry-ice. Use an aluminium foil to prevent loss of fragments during the cutting
 - or by lightly smashing the fragment in a mortar cooled by liquid nitrogen
 - 3) Put the fragment in one of the Potter homogenizers
 - 4) Add 9 volumes of « Mannitol buffer » (270 μL for 30 mg fragment for example)
 - 5) Homogenize at $0-6^{\circ}\text{C}$ (cold room or ice) with a glass pestle motor-driven at 500 to 1000 RPM. Stop as soon as the solution appears homogeneous, avoid over-homogenizing with subsequent heating of the solution. 3 passes are often sufficient.
 - 6) With a Pasteur pipette, transfer the solution into the 1.5 mL Eppendorf tube with the same number in the first series of Eppendorf tubes.
 - 7) Centrifuge at 650 g and $0-4^{\circ}\text{C}$ during 20 min
 - 8) **During centrifugation, put back to -80°C any remaining tissue fragment.**
 - 9) Transfer supernatant into the corresponding 1.5 mL Eppendorf tube in the second series of Eppendorf tubes.
 - 10) Resuspend the pellet with the volume of « Mannitol buffer » initially used, transfer the solution into the Potter homogenizer, and redo the previous steps of homogenization and centrifugation
 - 11) Pool the supernatant of the second centrifugation with the supernatant from the first one to obtain a 5% homogenate in "Mannitol buffer".
 - 12) The 5% homogenates will be kept at $0-6^{\circ}\text{C}$ for the assays.
- **Preparation of beef heart mitochondria analyzed in parallel to each series**
Mitochondrial beads have been provided to all centers of the French network by Prof. Joël Lunardi (CHU de Grenoble).
One bead is put in an Eppendorf tube and put at $0-4^{\circ}\text{C}$.
- When the mitochondrial suspension is melt, 10 μL of that suspension are transferred into 390 μL of "Mannitol buffer" (dilution 1/40), then 200 μL of the 1/40 solution are transferred to 800 μL of "Mannitol buffer". The dilution 1/200 is used for all the assays with the exception of complex III and complex IV assay for which a 1/800 dilution is used (50 μL of 1/200 dilution into 150 μL of "Mannitol buffer").
- **Preparation of cell suspension in "Mannitol buffer" (frozen pellets from cultured skin fibroblasts, or peripheral blood mononuclear cells)**
Resuspend the frozen pellet in "Mannitol buffer" (for example cells from a T75 flask are suspended in 200 μL of "Mannitol buffer")
Sonication during 5 seconds at moderate power is performed essentially to obtain an homogeneous solution.

- **Protein determination**

Prepare a 1/5 dilution in PBS of homogenates and cells suspension.

Measure the protein content of 10 μL of these 1/5 dilutions and of 10 μL of the 1/40 dilution of beef heart mitochondrial preparation using Pierce BCA kit and a 2 to 0.0625 mg/mL BSA standard curve.

Using "Mannitol buffer", adjust the protein concentration to 2 mg/mL for the homogenates and cells suspensions.

The protein concentration of the 1/40 dilution of the mitochondrial bead should be 1 mg/mL. If needed, modify the calculations of the rate of reaction.